



**SINO
MEDIC**

**N95 /N99 medical mask
Disposable medical masks**

disposable medical masks

- 
- **was established in Tainan, Taiwan in 2008**
 - **The company has long been engaged in the professional manufacturing of various medical masks for up to 12 years.**
 - **The marketing channels of medical masks in Taiwan: National Taiwan University Hospital, National Cheng Kung University Hospital, Kaohsiung Medical College, Large domestic clinics and franchised pharmacy across the country...etc.**
 - **During the anti-epidemic period of the COVID-19 in 2020, our factory was called by the government to become one of the team members of the 『 Taiwan National Mask Production Team 』 .**
 - **In 2020, the company extended its sales area to the international medical mask market. In the same year, the company obtained CE EU certification and FDA approval by the US FDA.**

Comparison of features

features Item	Our N95 (FFP2)	Well know branded N95
Filter method	The world's only "nanoporous filter" filter material that uses physical filtration.	1. Traditional electrostatic adsorption filter.
Pro.	<ol style="list-style-type: none"> 1. the PFE>99% is maintained for three years due to physical filtration. 2. Light and breath freely. 3. Washed more than 10 times , it still maintains BFE & PFE > 99 %. 4. Not afraid of moisture. 5. CE & FDA approval 	<ol style="list-style-type: none"> 1. The company is famous. 2. Head loop is comfortable to wear. 3. There are foundry factories all over the world. 4. CE & FDA(510 K) approval. 5. Widely known by doctors.
Con.	<ol style="list-style-type: none"> 1. Insufficient brand awareness internationally. 2. 6 million pieces per month 3. Niosh testing is accomplished, but waiting for certification 	<ol style="list-style-type: none"> 1. Afraid of moisture. 2. From the factory to the delivery and preservation process, if the electrostatic layer station ingon power is insufficient, will lead to a rapid decline in protection efficiency. 3. Only use once.
Price	<ol style="list-style-type: none"> 1. Stable price 2. High C/P value 3. Provide traders with high profits. 	<ol style="list-style-type: none"> 1. Floating price 2. High unit price 3. There are many fakes in the market

	Nano Technology N95 (FFP2)Medical Mask	Nano Technology N99 (FFP3)Medical Mask
Product number	<u>CMK-16</u>	<u>CMK-161</u>
washable	Can be washed	Can be washed
product photo		
specification	Mask size:11x15(±0.5CM) Nose piece length: 8(±0.5CM) Ear strap length: 22 (±0.5CM) Mask thickness: 0.62MM	Mask size:11x15(±0.5CM) Nose piece length: 8(±0.5CM) Ear strap length: 22 (±0.5CM) Mask thickness: 0.62MM
Quantity (Single Pack)	500pcs/carton	500pcs/carton
Carton size	46x53.5x39(CM)	46x53.5x39(CM)
Net weight of outer box	5.4KG	5.4KG
Gross weight of outer box	6.5KG	6.5KG
Manufacturing Origin	Taiwan(MIT)	Taiwan(MIT)
Shelf life	over 3 years	over 3 years
HS CODE	63079050115	
CE EU certification	✓	✓
FDA approved	✓	✓



**N95 (FFP2)
outward appearance**



**Internal appearance (With
sponge for nose pad)**



**S hook to enhanced
air tightness**



**Thickness less than
0.62 mm, ultra light
and breathable**



**Actual wearing
mask photo**

Are N95 (FFP2) medical masks there any other features of Nano Technology ?

- **Unique nano-filter material:** using the most top-quality materials, the high-efficiency filter material with many micro-holes, and the design with three-dimensional air chamber can achieve less burdensome air permeability, and is lighter and thinner than other brands, and breathe better.
- **Repeatable washing:** According to the "Taiwan Textile Research Institute (TTRI) Washing Test Report", after washing more than 10 times, the protection efficiency can still be maintained above PFE 99%, and the protection will not decrease after washing.
- **More economical and more environmentally friendly :** Take advantage of reproducible and washable, while prolonging the life of the mask, and also contribute to environmental protection.
- **Think about it:** how many people in the world throw away a piece of mask a day, and the inner meltblown layer is made of plastic material, which will never be decomposed all year round... The longer the time, the more pollution burden is added to the earth.

How to wash and clean the mask?

Cleaning Procedure:


1. It can be washed about 10 times.
2. After applying cleaning products with sterilization effects on both the inside and outside of the mask, use a toothbrush to gently clean the surface of the mask, rinse the mask directly under the tap with water on both sides, and then place it directly in a ventilated place to dry.

Caution:

- Remember not to twist or knead the body of the mask. This action will reduce the arrangement of the inner layer of the nanofilter material and affect the filtration efficiency.
- When you breathe, you smell the mask smelly /stinky or the surface of the mask is dirty / damaged. Please discard the mask immediately and replace it with a new one.

What is the difference between the NANO technical nano N95 (FFP2) medical mask and the medical masks on the market?

- General plane medical masks protect against particles PM3 that are larger than bacteria, but those harmful to the human body are all particles below PM2.5.
- Nanotechnology N95 (FFP2) medical mask can provide 99% protection against PM0.075 particles, providing more comprehensive protection.



1 Q & A for users

Is the unit price of NANO N95 mask more expensive than the KN95 made in China?

- **NANC N95 (FFP2) mask is the world's only and first new medical mask.**
- **It has stable protection performance under long-term storage, and has the air permeability when worn for a long time**
- **Tested by experiments, medical masks with traditional electrostatic filter media have been stored for more than 3 months, and the filtering power will drop by 20%~30%!**
- **Therefore, I hope that consumers can think seriously: Who wants to buy a mask that has only 80% filtration capacity remaining for 3 months?**
- **Or do you want to buy a medical mask that can maintain a stable 99% protection and more breathable and washable?**
- **A newly released new technology cannot be compared with the price of ordinary technology products on the market. What's more, the safety of life cannot be considered by price?**

2 Q & A for users

What are the certifications obtained for NANO N95 (FFP2) masks?

- **The factory has successfully obtained FDA approval with listing number & EUA and CE certification (B & C2).**
- **The International Bureau of Standards ISO9001 and ISO13485 certificates and orthographic marks.**
- **Obtained the only medical mask license in Taiwan that can manufacture N99 (FFP3).**
- **MIT smile badge expression [Safe, healthy and trustworthy] [Made in Taiwan, genuine products] A high standard certification mark.**

layer structure

The layer structure of Formosa Energy-Carbon FFP2 mask are as follows:

1. First layer: Waterproof material

PP water-repellent material is non-woven fabric, which can block the adhesion of pathogenic droplets or blood.

2. Second layer: Supportive fabric material

Special non-woven fabric support material can make the mask three-dimensional without deformation.

3. Third layer: Isolation fabric material


The PP isolation layer is non-woven fabric and can be blocked.

4. Fourth layer: High-tech nano-textile fabric material

High-tech nano-textile fabric filter layer can block bacteria and viruses by 99%.

5. Fifth layer: EP Bicomponent fabric material

Soft EP composite fiber non-woven fabric



A: From outside to inside, the specifications of the cloth are:

- 1: Waterproof layer: non-woven fabric
- 2: Supportive layer: PP non-woven fabric
- 3: Isolation layer: PP non-woven fabric
- 4: High-tech. filter layer: Nano-textile filter fabric
- 5: Soft composite layer: EP Soft fiber non-woven fabric

B: Nose bridge: Hot melt adhesive aluminum strip, 80x5

C: Ear strap: Polyester brocade spandex ribbon, 220x4

Strict quality control

- **Prevent infection between patients and medical staff.**
- **The elasticity of the earhook elastic band is good, and it will not cause ear pain after long wearing**
- **Both raw materials and production are made in Taiwan, and the production environment is carried out in a strictly controlled low-dust room.**
- **It has the characteristics of safety, hygiene, comfort and high protection.**

Taiwan Ministry of Health and Welfare Medical Equipment License



衛生福利部第一等級醫療器材許可證

衛部醫器製壹字第 0 [redacted] 號

中文名稱： [redacted] N95 醫用口罩(未滅菌)

英文名稱： [redacted]

類別： 第 I 類：一般及整型外科手術裝置 藥商名稱： [redacted]

規格： 空白 製造廠名稱： [redacted]
製造廠地址： [redacted]

效能： 限醫療器材管理辦法「醫療用衣物(I.4040)」第一等級鑑別範圍。

處方： 空白

前項醫療器材經本部審核與藥事法之規定相符應發給許可證以資證明

衛生福利部部長

陳時中



發證日期 107 年 12 月 12 日

有效日期 112 年 12 月 12 日

核准 展延 至	年 月 日	年 月 日	年 月 日	年 月 日
文號				

Taiwan Ministry of Health and Welfare Medical Equipment License



衛生福利部第一等級醫療器材許可證

衛部醫器製壹字第 0068 號

中文名稱： N99 醫用口罩 (未滅菌)

英文名稱：

類別： 第 I 類：一般及整型外科手術裝置 藥商名稱：

規格： 空白 製造廠名稱：
製造廠地址：

效能： 限醫療器材管理辦法「醫療用衣物(I.4040)」第一等級鑑別範圍。
處方： 空白

前項醫療器材經本部審核與藥事法之規定相符應發給許可證以資證明

衛生福利部部长

陳時中



發證日期 106 年 07 月 21 日

有效日期 111 年 07 月 21 日

核准 展延 至	年	月	日	年	月	日	年	月	日	年	月	日
	文號											

FDA(MSH)

CERTIFICATE OF **FDA** REGISTRATION 2020

This certifies that:

[Redacted]

is registered and has listed the following medical device with the U.S. Food and Drug Administration pursuant to Title 21, 807 et seq. of the United States Code of Federal Regulations:

Status: **Active**
Establishment Registration: [Redacted]
Owner/Operator Number: [Redacted]
Proprietary Name: [Redacted]
Device Classification Name: **RESPIRATOR, SURGICAL**
Product Code: **MSH**
Regulation Number: **878.4040**
Listing Number: [Redacted]

Compliance Engineering, Inc. makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. Compliance Engineering, Inc. assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."

The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. Compliance Engineering, Inc. is not affiliated with The U.S. Food and Drug Administration.

US Agent:
Compliance Engineering, Inc.
USA TEL: 516-984-7287
Taiwan TEL: 04-25685848

Kenneth Grant

Kenneth Grant
Executive Director
Compliance Engineering, Inc.

CE EN149 B+C2 certificate(1/8)



EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-P1

Respiratory protective devices, filtering half masks to protect against particles manufactured by

are tested and evaluated according to

EN 149:2001 + A1:2009 Respiratory Protective Devices - Filtering Half Masks to Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition

Brand Name: [REDACTED] Model: CMK-161
Filtering half mask
Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Ongoing successful performance in fulfilment of the requirements set out in **Personal Protective Equipment Regulation (EU) 2016/425** and harmonised standards, ensured by assessments based on **Annex 7 (Module C2) or Annex 8 (Module D)** of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on **30/06/2020** and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.




Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

UNIVERSAL

Verify the validity with the QR code



CE EN149 B+C2 certificate(2/8)



NB 2163

CERTIFICATE OF CONFORMANCE

Certificate No: 2163-

Respiratory protective devices, filtering half masks to protect against particles manufactured by

Continues to fulfil the requirements of

EN 149:2001 + A1:2009 Respiratory Protective Devices - Filtering Half Masks to Protect Against Particles - Requirements, Testing, Marking

Based on the evaluation of test reports and internal quality control audit reports according to EN 149+A1:2009 and Personal Protective Equipment Regulation (EU) 2016/425 Annex VII (Module C2). This certificate implies that the manufactured products show below are in conformance with the approved EU Type Examination model and meets the requirements of the regulation.

Product Definition


Model	Class	EU Type Examination Certificate		
		Serial No	Date	Issuing NB No
CMK-16	FFP2	216	30.06.2020	2163

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9**.
- Taking all measures necessary so that the manufacturing process and its monitoring ensure the homogeneity of production and conformity of the manufactured PPE with the type described in the EU type examination certificate.

This certificate is issued on **30/06/2020** and will be valid for one year, until **29/06/2021** if the manufacturer makes no major change in the product designs and manufacturing processes affecting the product performance on the essential health and safety requirement.




Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

Verify the validity with the QR code



CE EN149 B+C2 certificate(3/8)



TECHNICAL ASSESSMENT REPORT

REPORT DATE / NO: 30.06.2020 / 21

Manufacturer:
Address: No. 1:

This report is for the, given above, manufacturer prepared according to the test results obtained from Jiangsu Quality Supervision and Inspection Center for Special Safety Protection Products accredited by CNAS (China National Accreditation Service), signatory to ILAC MRA, with number L-7901 for the product identified below, dated 03.06.2020 with Serial Id STFWT202010821G based on EN 149: 2001 + A1: 2009 standard and the technical file dated 15 June 2020 Version 0 provided by the manufacturer. The sampling of the product is conducted under our supervision for testing from the manufacturing site of the client.

The technical file of the manufacturer, and risk evaluation against the essential health safety requirements and the test report evaluated for their relation with Essential Requirements of Personal Protective Equipment Regulation and found to be appropriate.

This report is an annex and an integral part of the EU Type Examination Certificate issued to the manufacturer. The test results and issued certificate belongs only to the tested model. The technical report consists of a total of 6 pages.

Product Description: Particle Filtering Half Mask
Classification: FFP2 NR
Trademark:



CE EN149 B+C2 certificate(4/8)



ESSENTIAL HEALTH and SAFETY REQUIREMENTS GIVEN IN EUROPEAN UNION REGULATION EU 2016/425 CORRESPONDING RISKS FOR THE PRODUCT

1.1. Design principles

1.1.1. Ergonomics

PPE must be so designed and manufactured that in the foreseeable conditions of use for which it is intended the user can perform the risk related activity normally whilst enjoying appropriate protection of the highest possible level. The test results with human subjects did not report any problem with the ergonomics of the product.

1.1.2. Levels and classes of protection

1.1.2.1. Highest level of protection possible

The optimum level of protection to be taken into account in the design is that beyond which the constraints by the wearing of the PPE would prevent its effective use during the period of exposure to the risk or normal performance of the activity.

1.1.2.2. Classes of protection appropriate to different levels of risk

Where differing foreseeable conditions of use are such that several levels of the same risk can be distinguished, appropriate classes of protection must be taken into account in the design of the PPE.

1.2. Innocuousness of PPE

1.2.1. Absence of risks and other inherent nuisance factors

PPE must be so designed and manufactured as to preclude risks and other nuisance factors under foreseeable conditions of use. The manufacturer declares in his technical file that according to the results of risk analysis and the material properties they use in the manufacturing, the product has no hazardous content for health.

1.2.1.1. Suitable constituent materials

The materials of which the PPE is made, including any of their possible decomposition products, must not adversely affect the health or safety of users. The material selection is processed in the technical manufacturing process and documented.

1.2.1.2. Satisfactory surface condition of all PPE parts in contact with the user

Any part of the PPE that is in contact or is liable to come into contact with the user when the PPE is worn must be free of rough surfaces, sharp edges, sharp points and the like which could cause excessive irritation or injuries is evaluated and reported in the test report.

1.2.1.3. Maximum permissible user impediment

Any impediment caused by PPE to movements to be made, postures to be adopted and sensory perception must be minimized; nor must PPE cause movements which endanger the user or other persons.

1.3. Comfort and effectiveness

1.3.1. Adaptation of PPE to user morphology

PPE must be designed and manufactured in such a way as to facilitate its correct positioning on the user and to remain in place for the foreseeable period of use, bearing in mind ambient factors, the actions to be carried out and the postures to be adopted. For this purpose, it must be possible to adapt the PPE to fit the morphology of the user by all appropriate means, such as adequate adjustment and attachment systems or the provision of an adequate range of sizes.

1.3.2. Lightness and design strength

PPE must be as light as possible without prejudicing design strength and efficiency.

Apart from the specific additional requirements which they must satisfy in order to provide adequate protection against the risks in question (see 3), PPE must be capable of withstanding the effects of ambient phenomena inherent under the foreseeable conditions of use

1.4. Information supplied by the manufacturer

The notes that must be drawn up by the former and supplied when PPE is placed on the market must contain all relevant information on:

- In addition to the name and address of the manufacturer and/or his authorized representative established in the Community
- Storage, use, cleaning, maintenance, servicing and disinfection, cleaning, maintenance or disinfectant protection recommended by manufacturers must have no adverse effect on PPE or users when applied in accordance with the relevant instructions;
- Performance as recorded during technical tests to check the levels or classes of protection provided by the PPE in question;
- Suitable PPE accessories and the characteristics of appropriate spare parts;
- The classes of protection appropriate to different levels of risk and the corresponding limits of use;
- The obsolescence deadline/period of obsolescence of PPE or certain of its components;
- The type of packaging suitable for transport;
- The significance of any markings (see 2.12)
- Where appropriate the references of the Directives applied in accordance with Article 5(6) (b);
- The name, address and identification number of the notified body involved in the design stage of the PPE

These notes, which must be precise and comprehensible, must be provided at least in the official language(s) of the member state of destination



CE EN149 B+C2 certificate(5/8)



2. ADDITIONAL REQUIREMENTS COMMON TO SEVERAL CLASSES OR TYPES OF PPE

2.1. PPE incorporating adjustment systems

If PPE incorporates adjustment systems, the latter must be designed and manufactured so that, after adjustment, they do not become undone unintentionally in the foreseeable conditions of use.

2.3. PPE for the face, eyes and respiratory system

Any restriction of the user's face, eyes, field of vision or respiratory system by the PPE shall be minimised.

The screens for those types of PPE must have a degree of optical neutrality that is compatible with the degree of precision and the duration of the activities of the user.

If necessary, such PPE must be treated or provided with means to prevent misting-up.

Models of PPE intended for users requiring sight correction must be compatible with the wearing of spectacles or contact lenses.

2.4. PPE subject to ageing

If it is known that the design performance of new PPE may be significantly affected by ageing, the month and year of manufacture and/or, if possible, the month and year of obsolescence must be indelibly and unambiguously marked on each item of PPE placed on the market and on its packaging.

If the manufacturer is unable to give an undertaking with regard to the useful life of the PPE, his instructions must provide all the information necessary to enable the purchaser or user to establish a reasonable obsolescence month and year, taking into account the quality level of the model and the effective conditions of storage, use, cleaning, servicing and maintenance.

Where appreciable and rapid deterioration in PPE performance is likely to be caused by ageing resulting from the periodic use of a cleaning process recommended by the manufacturer, the latter must, if possible, affix a marking to each item of PPE placed on the market indicating the maximum number of cleaning operations that may be carried out before the equipment needs to be inspected or discarded. Where such a marking is not affixed, the manufacturer must give that information in his instructions. The product is for single use and tested with simulated wearing conditioning.

2.6. PPE for use in potentially explosive atmospheres

PPE intended for use in potentially explosive atmospheres must be designed and manufactured in such a way that it cannot be the source of an electric, electrostatic or impact-induced arc or spark likely to cause an explosive mixture to ignite.

2.8. PPE for intervention in very dangerous situations

The instructions supplied by the manufacturer with PPE for intervention in very dangerous situations must include, in particular, data intended for competent, trained persons who are qualified to interpret them and ensure their application by the user.

The instructions must also describe the procedure to be adopted in order to verify that PPE is correctly adjusted and functional when worn by the user.

Where PPE incorporates an alarm which is activated in the absence of the level of protection normally provided, the alarm must be designed and placed so that it can be perceived by the user in the foreseeable conditions of use.

2.9. PPE incorporating components which can be adjusted or removed by the user

Where PPE incorporates components which can be attached, adjusted or removed by the user for replacement purposes, such components must be designed and manufactured so that they can be easily attached, adjusted and removed without tools.

2.12. PPE bearing one or more identification or recognition marks directly or indirectly relating to health and safety

The identification or recognition marks directly or indirectly relating to health and safety affixed to these types or classes of must preferably take the form of harmonized pictograms or ideograms and must remain perfectly legible throughout the foreseeable useful life of the PPE. In addition, these marks must be complete, precise and comprehensible so as to prevent any misinterpretation; in particular, where such marks incorporate words or sentences, the latter must appear in the official language(s) of the Member State where the equipment is to be used.

If PPE (or a PPE component) is too small to allow all or part of the necessary marking to be affixed, the relevant information must be mentioned on the packing and in the manufacturer's notes.

3. ADDITIONAL REQUIREMENTS SPECIFIC TO PARTICULAR RISKS

3.10.1. Respiratory protection

PPE intended for the protection of the respiratory system must make it possible to supply the user with breathable air when exposed to a polluted atmosphere and/or an atmosphere having an inadequate oxygen concentration.

The breathable air supplied to the user by PPE must be obtained by appropriate means, for example after filtration of the polluted air through PPE or by supply from an external unpolluted source.

The constituent materials and other components of those types of PPE must be chosen or designed and incorporated so as to ensure appropriate user respiration and respiratory hygiene for the period of wear concerned under the foreseeable conditions of use.

The leak-tightness of the facepiece and the pressure drop on inspiration and, in the case of the filtering devices, purification capacity must keep contaminant penetration from a polluted atmosphere low enough not to be prejudicial to the health or hygiene of the user.

The PPE must bear details of the specific characteristics of the equipment which, in conjunction with the instructions, enable a trained and qualified user to employ the PPE correctly.

In the case of filtering equipment, the manufacturer's instructions must also indicate the time limit for the storage of new filters kept in their original packaging.



CE EN149 B+C2 certificate(6/8)



Technical Assessment of EN 149: 2001 + A1: 2009 Standard and other Standards it refers to. Clauses Corresponding to the (EU) 2016/425 Regulation. Essential Health and Safety Requirements given above.

Conforming to EN 149:2001 + A1:2009 Standard Requirements

Article 5	<p>Classification: Particle Filtering Half Mask</p> <p>The mask subject to evaluation based on the test results and technical file provided by the manufacturer is classified as: Filtering Efficiency and maximum Total Inward Leakage: Classified as FFP2 Mask is classified for single shift use, NR</p>																																					
Article 7.4	<p>Packing: Particle filtering half masks are packaged to protect them from contamination before use and with cardboard boxes to prevent mechanical damage, the masks are in plastic sealed bags in the card box. The packaging design and the product is considered to withstand the foreseeable conditions of use based on the visual inspection results given in the test report.</p>																																					
Article 7.5	<p>Material: Materials used in particle filtering half masks, according to the simulated wearing treatment and temperature conditioning results; It is understood it withstands handling and wear over the period for which the particle filtering half mask is designed to be used, it suffered mechanical failure of the facepiece or straps, any material from the filter media released by the air flow through the filter has not constitute a hazard or nuisance for the wearer. The manufacturer declares that the materials used in manufacturing of the mask does not have an adverse affect the health and safety of users.</p> <p>Based on the test results, the masks did not collapse when subject to simulated wearing and temarature conditioning. No nuisance situation is reported during the practical performance tests by human subjects.</p>																																					
Article 7.6	<p>Cleaning and Disinfection: Particle filtering half mask is not designed to be as re-usable. No cleaning or disinfection procedure provided by the manufacturer.</p>																																					
Article 7.7	<p>Practical Performance :</p> <p>The test report indicates that the human subjects did not face any difficulty in performing the exercises while they were weared by the sample masks, in walking test or work simulation tests. The wearers did not report any failure by means of head harness / straps/ earloops comfort, security of fastenings and field of vision. Also no imperfections reported during total inward tests about the comfort, field of vision and fastening issues.</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th>Assessed Elements</th> <th>Positive</th> <th>Negative</th> <th>Requirements in accordance with EN 149:2001 + A1:2009 and Result</th> </tr> </thead> <tbody> <tr> <td>2.Head harness comfort</td> <td>2</td> <td>0</td> <td rowspan="3">Positive results are obtained from the test subjects No imperfections</td> </tr> <tr> <td>3.Security of fastenings</td> <td>2</td> <td>0</td> </tr> <tr> <td>5.Field of vision</td> <td>2</td> <td>0</td> </tr> </tbody> </table> <p>Conditioning : (A.R.) As Received, original</p>	Assessed Elements	Positive	Negative	Requirements in accordance with EN 149:2001 + A1:2009 and Result	2.Head harness comfort	2	0	Positive results are obtained from the test subjects No imperfections	3.Security of fastenings	2	0	5.Field of vision	2	0																							
Assessed Elements	Positive	Negative	Requirements in accordance with EN 149:2001 + A1:2009 and Result																																			
2.Head harness comfort	2	0	Positive results are obtained from the test subjects No imperfections																																			
3.Security of fastenings	2	0																																				
5.Field of vision	2	0																																				
Article 7.8	<p>Finish of Parts: The test report states that the particle filtering half masks, which are likely to come into contact with the user, do not have sharp edges and do not contain burrs.</p>																																					
Article 7.9.1	<p>Total Inward Leakage:</p> <p>The Total Inward Leakage test is conducted by 10 individual in an aerosol chamber with a walking band, and samples are taken during the conduction of the exercises defined in the standard. The samples used in the test are subjected to the conditioning required in the standard as Temperature conditioning and as received. The face dimensions of the subjects are also reported. The measurement details for each subject and for each exercsize are available in the test report.</p> <p>It was reported that; At least 47 out of 50 exercise measurement results are smaller or equal to 11%. At least 9 of 10 individual's arithmetic mean is smaller or equal to 8%.</p> <p style="text-align: center;">According to the reported results, the product meets the limits for FFP1 and FFP2 classification.</p>																																					
Article 7.9.2	<p>Penetration of filter material: Sodium Chloride Testing</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th>Condition</th> <th>No. of Sample</th> <th>Sodium Chloride Testing 95 L/min max (%)</th> <th>Requirements in accordance with EN 149:2001 + A1:2009</th> <th>Result</th> </tr> </thead> <tbody> <tr> <td>(A.R.)</td> <td>19#</td> <td>0.68</td> <td rowspan="3">FFP1 ≤ 20 %</td> <td rowspan="9">Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 classes.</td> </tr> <tr> <td>(A.R.)</td> <td>20#</td> <td>0.69</td> </tr> <tr> <td>(A.R.)</td> <td>21#</td> <td>0.71</td> </tr> <tr> <td>(S.W.)</td> <td>22#</td> <td>0.89</td> <td rowspan="2">FFP2 ≤ 6 %</td> </tr> <tr> <td>(S.W.)</td> <td>23#</td> <td>0.94</td> </tr> <tr> <td>(S.W.)</td> <td>24#</td> <td>0.91</td> <td rowspan="3">FFP3 ≤ 1 %</td> </tr> <tr> <td>(M.S. T.C.)</td> <td>25#</td> <td>1.02</td> </tr> <tr> <td>(M.S. T.C.)</td> <td>26#</td> <td>1.11</td> </tr> <tr> <td>(M.S. T.C.)</td> <td>27#</td> <td>1.07</td> <td></td> </tr> </tbody> </table> <p>Conditioning : (M.S.) Mechanical Strength (T.C.) Temperature Conditioning (A.R.) As Received, original (S.W.) Simulated wearing treatment</p> <p style="text-align: right;">95 L/min = 1.6 dm³.sn⁻¹</p>	Condition	No. of Sample	Sodium Chloride Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result	(A.R.)	19#	0.68	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 classes.	(A.R.)	20#	0.69	(A.R.)	21#	0.71	(S.W.)	22#	0.89	FFP2 ≤ 6 %	(S.W.)	23#	0.94	(S.W.)	24#	0.91	FFP3 ≤ 1 %	(M.S. T.C.)	25#	1.02	(M.S. T.C.)	26#	1.11	(M.S. T.C.)	27#	1.07	
Condition	No. of Sample	Sodium Chloride Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result																																		
(A.R.)	19#	0.68	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 classes.																																		
(A.R.)	20#	0.69																																				
(A.R.)	21#	0.71																																				
(S.W.)	22#	0.89	FFP2 ≤ 6 %																																			
(S.W.)	23#	0.94																																				
(S.W.)	24#	0.91	FFP3 ≤ 1 %																																			
(M.S. T.C.)	25#	1.02																																				
(M.S. T.C.)	26#	1.11																																				
(M.S. T.C.)	27#	1.07																																				



CE EN149 B+C2 certificate(7/8)



Penetration of filter material : Paraffin Oil Testing						
Condition	No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result		
(A.R.)	28#	0.98	FFP1 ≤ 20 % FFP2 ≤ 6 % FFP3 ≤ 1 %	Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 classes.		
(A.R.)	29#	1.05				
(A.R.)	30#	1.03				
(S.W.)	31#	1.25				
(S.W.)	32#	1.19				
(S.W.)	33#	1.21				
(M.S. T.C.)	34#	1.33				
(M.S. T.C.)	35#	1.46				
(M.S. T.C.)	36#	1.39				
Conditioning : (M.S.) Mechanical Strength (T.C.) Temperature Conditioning (A.R.) As Received, original (S.W.) Simulated wearing treatment						
Article 7.10	Compatibility with skin: In Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other adverse effect on health was not reported. (No negative reporting on practical performance and TIL test results)					
Flammability :						
Article 7.11	Condition	No. of Sample	Visual inspection	Requirements in accordance with EN 149:2001 + A1:2009	Result	
	(A.R.)	37#	Didn't burn	Filtering half mask shall not burn or not continue to burn for more than 5 s after removal from the flame	Passed Laboratory claims that the tested items did not burn for 5 seconds and fulfils the requirement of the standard	
	(A.R.)	38#	Didn't burn			
	(T.C.)	39#	Didn't burn			
	(T.C.)	40#	Didn't burn			
Conditioning : (A.R.) As Received, original (T.C.) Temperature Conditioning						
Carbon dioxide content of the inhalation air:						
Article 7.12	Condition	No. of Sample	CO ₂ content of the inhalation air [%] by volume	An average CO ₂ content of the inhalation air	Requirements in accordance with EN 149:2001 + A1:2009	Result
	(A.R.)	41#	0.55	0,56 [%]	CO ₂ content of the inhalation air shall not exceed an average of 1,0% by volume	Passed Filtering half masks fulfil requirements of the standard
	(A.R.)	42#	0.57			
	(A.R.)	43#	0.56			
Conditioning : (A.R.) As Received, original						
Article 7.13	Head harness: In Practical Performance and TIL test reports no adverse effects have been reported for donning and remove of the mask also the results of these tests indicates that the head harness are capable of holding the mask firmly enough.					
Article 7.14	Field of vision: In Practical Performance report, no adverse effects were reported for the field of vision availability when the mask is worn.					
Article 7.15	Exhalation Valve(s): The model under inspection have no valves.					
Article 7.16	Breathing Resistance: Inhalation The overall evaluation of the results gathered for 9 different samples 3 as received, 3 with temperature conditioning, 3 simulated wearing treatment complies with the limits given in the standard for FFP1, FFP2 and FFP3 classes. This is valid for inhalation results for 30 L/min, 95 L/min and exhalation at 160 L/min. The measurement details for each single mask tested are available in the test report. Passed.					

CE EN149 B+C2 certificate(8/8)



Article 7.17	Clogging: This test is not applied to Particle Filtering Half Mask which is not reusable. <i>(For single shift use devices, the clogging test is optional test. For re-usable devices test is mandatory.)</i>
Article 7.18	Demountable Parts: There are no demountable parts of the mask.
Article 8	Testing: All tests conducted according to Clause 8 of this standard is available in the test report and are evaluated in this report for qualification and classification of the mask.
Article 9	Marking – Packaging: Necessary markings are available on the product package (box). The manufacturer and its trademark is clearly visible. The type of the mask and the classification including the status of re-usability, the reference to EN 149:2001+A1:2009 standard, the end date of shelf life, using and storage instructions and pictograms and CE mark are available on the product package. The above evaluation is based on the technical document for packaging and marking, for box design. Verified on the Annex 5 of the technical file. The technical documentation for mask design (drawing) also evaluated for marking requirements, drawing CMK-1619. The mask template (drawing) indicates that the mask will carry information about the manufacturer / trademark (Formosa Energy-Carbon) of the manufacturer. Type of mask, the reference to EN 149+A1:2009 standard and classification including the re-usability of the mask. The manufacturer also printed CE mark with our Notified Body number. The mask do not have sub-assemblies. Even the tested sample by the laboratory do not carry necessary marking information as stated in the technical documentation, the manufacturer shall follow marking instructions for serial production. Model drawing CMK-1619 exists in the technical file of the manufacturer, Annex 4 of technical file.
Article 10	Information to be supplied by the manufacturer: In each of the smallest commercially available packaging of the product, implementation (installation instructions) pre-use controls, warning and usage limitations, storage and meanings of symbols / pictograms are defined. User instruction document in the technical file found to be appropriate. The manufacturer shall include this documented user information text in every smallest commercially available package, Annex 2 of Technical file.

PREPARED BY	APPROVED BY
Neslihan EKE BİRTÜRK PPE Expert	Suat KAÇMAZ General Manager



ISO9001 & ISO13485 Certificate



THE INTERNATIONAL CERTIFICATION NETWORK

CERTIFICATE

CISQ/IMQ as an IQNet Partner hereby states that the organization

for the following scope:

Design and manufacture of disposable respirators

Further clarifications regarding the applicability of ISO 9001:2015 requirements may be obtained by consulting the organization

has implemented and maintains a
Quality Management System

which fulfills the requirements of the following standard

ISO 9001:2015

Issued on: 2017-08-03

First issued on: 2017-08-03

for the validity date, please refer to the original certificate* issued by IMQ

Registration Number: IT - 112095



Alex Stoichitoui
President of IQNET



Ing. Claudio Provetti
President of CISQ

IQNet Partners:**

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil
FONDONORMA Venezuela ICONTEC Colombia Inspecta Certification Finland INTECO Costa Rica
IRAM Argentina IQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland PCBC Poland
Quality Austria Austria RR Russia SIGE México SII Israel SIQ Slovenia SIRIM QAS International Malaysia
SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey Vinçotte Belgium YUQS Serbia
IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

* This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document

** The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under www.iqnet-certification.com



THE INTERNATIONAL CERTIFICATION NETWORK

CERTIFICATE

CISQ/IMQ as an IQNet Partner hereby states that the organization

for the following scope:

Design and manufacture of disposable respirators

Further clarifications regarding the applicability of ISO 13485:2016 requirements may be obtained by consulting the organization

has implemented and maintains a
Quality Management System

which fulfills the requirements of the following standard

ISO 13485:2016

Issued on: 2017-08-03

First issued on: 2017-08-03

for the validity date, please refer to the original certificate* issued by IMQ

Registration Number: IT - 112127



Alex Stoichitoui
President of IQNET



Ing. Claudio Provetti
President of CISQ

IQNet Partners:**

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil
FONDONORMA Venezuela ICONTEC Colombia Inspecta Certification Finland INTECO Costa Rica
IRAM Argentina IQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland PCBC Poland
Quality Austria Austria RR Russia SIGE México SII Israel SIQ Slovenia SIRIM QAS International Malaysia
SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey Vinçotte Belgium YUQS Serbia
IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

* This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document

** The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under www.iqnet-certification.com

「0.075 μ m filtration rate」 Test Report

劉建新 (06-243-6135)

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



TEST REPORT TUCHENG

Date: May.07,2020 Date of Receipt: Apr.30,2020

Report No.: TFF9D786 Quantity: 1PC Page Order/Pages: (P1/3) Ref. No.: NIL

Report Title: Item: Mask

Address:

Test Items		Test Results	Test Methods
Protection Efficiency of Mask(%) (0.075 μ m NaCl,CMD)	1	98.72	CNS 14755 Z2125-2011 Flow rate:85.1 (liter/min) <i>This test report is for reference only and should not be used as a formal document.</i>
	2	99.62	
	3	98.87	
	4	98.95	
	5	98.69	
	6	98.88	
	7	99.47	
	8	99.04	
	9	99.27	
	10	99.05	
	Ave.	99.06	

Test Results
98.69%~
99.62%

Note: Sample description was given by the client: TWEC N95 Respirator
LOT: 20200428

Note: As the remaining sample was asked to return along with the test report, re-testing would not be possible.

Note: The test report is the test results issued by the testing institution as requested by the consignor, it shall not determine the legitimacy of the product.

Protection Efficiency of Mask(%)

- Note: 1.This report is only responsible for the submitted sample(s), which will be kept for one month period.
2.This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.
3.The test report should not be used for public advertisement and commercial promotion.

Authorized by president of
Taiwan Textile Research Institute

Jui-hung kao

Director,
Department of Testing and Certification

Department of Testing and Certification, Taiwan Textile Research Institute
No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R.O.C.)
Tel : +886-2-22670321 ext. 7107, 7110
Fax : +886-2-22675108 , +886-2-22689839

「Respiratory Impedance」 Test Report

劉建新 (06-243-6135)

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



TEST REPORT TUCHENG

Date: May.07,2020 Date of Receipt: Apr.30,2020
 Report No.: TFF9D786 Quantity: 1PC Page Order/Pages: (P2/3) Ref. No.: NIL
 Report Title: _____ Item: Mask
 Address: _____

Test Items		Test Results	Test Methods
Inhalation Resistance (Pa)	1	101.92	CNS 14755 Z2125-2011 Flow rate:85.1 (L/min)
	2	97.02	
	3	89.18	
	4	95.06	
	5	93.10	
	6	101.92	
	8	104.86	
	10	99.96	
	Ave.	102.90	
	10	92.12	
	10	97.80	
Exhalation Resistance (Pa)	1	100.94	CNS 14755 Z2125-2011 Flow rate:85.1 (L/min)
	2	85.26	
	3	85.26	
	4	87.22	
	5	79.38	
	6	84.28	
	8	97.02	
	10	91.14	
	10	92.12	
	10	86.24	
	Ave.	88.89	

Inhalation Resistance

Exhalation Resistance

Test Results: Higher than standard value

This test report is for reference only and should not be used as a formal document.

Note: $ImmH2O=9.8Pa$.

Note: Sample description was given by the client: TWEC N95 Respirator

LOT: 20200428

Note: As the remaining sample was asked to return along with the test report, re-testing would not be possible.

Note: The test report is the test results issued by the testing institution as requested by the consignor, it shall not determine the legitimacy of the product.

- Note: 1.This report is only responsible for the submitted sample(s), which will be kept for one month period.
 2.This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.
 3.The test report should not be used for public advertisement and commercial promotion.

Authorized by president of **Jui-hung Kao**
 Taiwan Textile Research Institute
 Director, Department of Testing and Certification

Department of Testing and Certification, Taiwan Textile Research Institute
 No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R.O.C.)
 Tel : +886-2-22670321 ext. 7107, 7110
 Fax : +886-2-22675108 , +886-2-22689839

「BFE」 bacterial filtration rate

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



1

正 本
ORIGINAL

試驗報告 土城場區
TEST REPORT TUCHENG

日期 Date: May. 22, 2017 收件日期 Date of Receipt: May. 01, 2017
 報告編號 Report No.: TFF6E003-D 數量 Quantity: IPC 報告頁次/頁數 Page Order/Pages: (P1/1) 來文字號 Ref. No.: NIL
 報告抬頭 Report Title: _____ 試件類別 Item: MASK
 地址 Address: NO. _____

Test Items	Test Results	Test Methods
Bacterial Filtration	1 > 99.9	CNS 14774 T5017-2011 9.2
Efficiency (BFE)(%)	2 > 99.9	CNS 14775 T4037-2003
Staphylococcus aureus	3 > 99.9	
ATCC 6538	4 > 99.9	
	5 > 99.9	

Test Results:
99.9%

Note: Control average: 2016 CFU.
 Note: Mean particle size: 3.0 μm.
 Note: Testing side: inside of specimen.
 Note: Testing area: 39.5 cm².
 Note: As requested by the client, test location: A flat area of mask.
 Note: Sample description is given by the client: FORMOSA ENERGY-CARBON N95 RESPIRATOR
 Note: As the remaining sample was asked to return along with the test report, re-testing would not be possible.
 Note: The test report is the test results issued by the testing institution as requested by the consignor, it shall not determine the legitimacy of the product.

Bacterial Filtration Efficiency(BFE)

附記: 1.本報告僅對樣品負責, 樣品保留期限為一個月。
 Note: This report is only responsible for the submitted sample(s), which will be kept for one month period.
 2.本報告非經本檢測及驗證部書面同意, 不得摘錄複製, 但全部複製除外。
 This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.
 3.試驗報告所載事項不得作為公開廣告及商業推銷之用。
 The test report should not be used for public advertisement and commercial promotion.



Guo-Tsun Jou
Director.

Authorized by president
 Taiwan Textile Research Institute
Department of Testing and Certification
 紡織所檢測及驗證部 新北市土城區承天路六號
 Department of Testing and Certification, Taiwan Textile Research Institute
 No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
 Tel : +886-2-22670321 ext. 7107, 7110
 +886-2-22698229, +886-2-22698230
 Fax: +886-2-22675108, +886-2-22689839

「Synthetic Blood Penetration」 Test Report

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



正 本
ORIGINAL

試驗報告 土城場區
TEST REPORT TUCHENG

日期: May. 22, 2017 收件日期: May. 01, 2017
 報告編號: TFF6E003-C 數量: IPC 報告頁次/頁數: (P1/1) 來文字號: NIL
 報告抬頭: [Redacted] 試件類別: MASK
 地址: [Redacted]

Test Items	Test Results	Test Methods
Synthetic Blood Penetration Pressure: 80 mmHg	1	none
	2	none
	3	none
	4	none
	5	none
	6	none
	7	none
	8	none
	9	none
	10	none

**Test Results:
Sputtering
pressure
through
blood 120**

**Synthetic Blood
Penetration
Pressure**

Note: Sample description is given by the client: FORMOSA ENERGY-CARBON N95 RESPIRATOR

Note: As the remaining sample was asked to return along with the test report, it should be tested as possible.

Note: The test report is the test results issued by the testing institution as requested by the consignor, it shall not determine the legitimacy of the product.

附記: 1.本報告僅對樣品負責, 樣品保留期限為一個月。
 Note: This report is only responsible for the submitted sample(s), which will be kept for one month period.
 2.本報告非經本檢測及驗證部書面同意, 不得摘錄複製, 但全部複製除外。
 This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.
 3.試驗報告所載事項不得作為公開廣告及商業推銷之用。
 The test report should not be used for public advertisement and commercial promotion.



Gwo-Tsun Jow

財團法人紡織產業綜合研究所
 所長授權核發人: Director.
 Authorized by president of Department of Testing and
 Taiwan Textile Research Institute Certification

紡織所檢測及驗證部 新北市土城區承天路六號
 Department of Testing and Certification, Taiwan Textile Research Institute
 No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
 Tel : +886-2-22670321 ext. 7107, 7110
 +886-2-22698229, +886-2-22698230
 Fax : +886-2-22675108, +886-2-22689839

「Repeat 10 Washes」 Test Report

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



日期 Date: 2018.12.26 收件日期 Date of Receipt: 2018.10.05 試驗報告 TEST REPORT 土城場區 TUCHENG 正本 ORIGINAL

報告編號 Report No.: TFF71109 數量 Quantity: 1件 報告頁次/頁數 Page Order/Pages: (P1/2) 來文字號 Ref. No.: 空白

報告抬頭 Report Title: 試件類別 Item: 口罩

地址 Address: _____

試驗項目	試驗結果	試驗方法
防護效率(%) (0.075 μm NaCl, CMD)	1	99.01
	2	99.24
	3	99.94
	4	99.77
	5	99.76
	6	99.79
	7	99.75
	8	99.85
	9	99.56
	10	99.47
	平均	99.61

Test Results:
99.01%~
99.94%

註：水洗方法：CNS 13752 L3243-1996 A法 篩乾x10

註：依委託者所提供來樣資料為：台灣精碳N95醫用口罩

註：試驗報告僅就委託者之委託事項提供試驗結果，不對產品合法性做判斷。

**Submicron particles
0.075 μm Filtration rate**



「Repeat 10 Washes」 Test Report

TTRI 財團法人紡織產業綜合研究所
Taiwan Textile Research Institute



日期 Date: 2018.12.26 收件日期 Date of Receipt: 2018.10.05 試驗報告 TEST REPORT 土城場區 TUCHENG 正本 ORIGINAL
 報告編號 Report No.: TFF7J109 數量 Quantity: 1件 報告頁次/頁數 Page Order/Pages: (P2/2) 來文字號 Ref. No.: 空 白
 報告抬頭 Report Title: 試件類別 Item: 口罩
 地址 Address:

試驗項目	試驗結果	試驗方法
吸氣阻抗(Pa)	1	149.94
	2	139.16
	3	154.84
	4	134.26
	5	133.28
	6	149.94
	7	122.50
	8	131.32
	9	137.20
	10	130.34
	平均	138.28
呼氣阻抗(Pa)	1	152.88
	2	148.96
	3	175.42
	4	149.94
	5	144.06
	6	168.56
	7	121.52
	8	147.00
	9	130.34
	10	151.90
	平均	149.06

Inhalation Resistance

Exhalation Resistance

Test Results: Higher than standard value

註：水洗方法：CNS 13752 L3243-1996 A/a 篩乾x10

註：1mmH2O=9.8Pa

註：依委託者所提供來樣資料為：台灣精碳N95醫用口罩

註：試驗報告僅就委託者之委託事項提供試驗結果，不對產品合法性做判斷。

附記：1.本報告僅對樣品負責，樣品保留期限為一月
 Note: This report is only responsible for the sample(s), which will be kept for one month
 2.本報告非經本檢測及驗證部書面同意，不得全部複製除外。
 This report cannot be reproduced in any full context, without the prior approval of this Department of Testing and Certification
 3.試驗報告所載事項不得作為公開廣告及商業用途
 The test report should not be used for public advertisement and commercial promotion.

財團法人紡織產業綜合研究所
 所長授權核發人：
 Authorized by president of
 Taiwan Textile Research Institute

 財團法人紡織產業綜合研究所檢測及驗證部 新北市土城區承天路6號
 Department of Testing and Certification, Taiwan Textile Research Institute
 No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R.O.C.)
 Tel : +886-2-22670321 ext. 7107, 7110
 Fax : +886-2-22675108, +886-2-22689839

檢測及驗證部 周國村

Front view of the Package (Chinese)

防塵や医療用で 国に認定された最高レベルのN95マスク

+ 國家雙認證高機能N95最高級口罩

N95

空汙 健康 流行病 運動 出差 工作

抗菌防護PM2.5

高科技奈米五層防護層

通過歐盟及台灣多項認證

防曬抗UV達99%紫外線

機車、通勤族 家庭主婦 油煙攤販 油漆工 保護肺部最佳的選擇

可水洗重複使用

花粉 20 μ m
PM2.5 2.5 μ m
水霧 5 μ m
酸鹼 10 μ m
細菌 0.5 μ m

石粉 0.09 μ m
油煙 1 μ m
研磨粉塵 1 μ m

0.075 μ m 防護等級

呼吸好空氣

145cm

超大氣流量
呼吸好舒暢

特殊防水設計
病毒防護等級再升級

高透氣、安全、衛生、舒適、最高防護、不悶熱

75mm
105mm
30mm

190mm

145mm

Back view of the Package (Chinese)

38mm
91mm
5mm

對人體有害的為PM2.5以下的細懸浮微粒

真正對人體有害的粒子

花粉 20µm	塵蟎 10µm	灰塵 5µm	細菌 3µm	PM2.5 2.5µm	油煙/香菸 1µm	麵粉 1µm	流感病毒 0.3µm	石棉 0.09µm	國家認證 標準 0.075µm

奈米科技N95口罩過濾效能 **PFE99%**

使用方法

-
-
-
-
-

1 取出口罩 2 將耳帶掛於兩耳後 3 下壓鼻片
4 調整口罩至最密位置 5 將5勾確實勾緊兩耳帶使口罩緊密貼合於面部

名稱：
適用對象：成人
規格：1. 立體耳掛式
2. 口罩尺寸：11.5cm(±0.5cm)X11.5cm(±0.5cm)
3. 耳帶長度：22cm(±0.5cm)
4. 鼻片長度：8cm(±0.5cm)
數量：1入(枚)
材料：PP不織布、特殊不織布支撐材、防護濾網、
高防護奈米科技濾網、複合超柔軟親膚層
批號：9520191101
製造日期：2019/11/01
保存期限：五年
許可字號：
製造商：
地址：
進口商/經銷商：

建議售價NTD150

4 719871 810866

190mm

145mm

Front view of package (US regulations)



190mm

145mm

Front view of package (US regulations)



190mm

145mm

Front view of packaging (European regulations)

The image shows the front view of a respirator packaging box. The box is primarily green and white. At the top right, there is a gold shield-shaped logo with the text "FFP2". Below this, there are five circular icons representing different environments: Agriculture, Healthcare, Home, Industrial, and Work. The main text on the box reads "The best option to protect the lungs" in large yellow letters. To the right of this text, there are three white boxes with black text: "99% UV Resistance", "Antibacterial Protection PM2.5", and "Nanotechnology 5 protective layers". Below these boxes are several certification logos: CE, TAP, IAT, ISO, and others. The bottom section of the box features a large illustration of a hand wearing a respirator. To the left of the hand, there is a circular diagram showing various particles being filtered out, including Pollen, Asbestos, Dust, and others. Below this diagram, there are two blue boxes with white text: "Breathable Fabric" and "Waterproof Design- Upgrading protection from viruses". At the bottom of the box, there is a red stamp that says "MADE IN TAIWAN". The overall design is clean and professional, emphasizing the product's safety and effectiveness.

FFP2

Agriculture Healthcare Home Industrial Work

The best option to protect the lungs

99% UV Resistance

Antibacterial Protection PM2.5

Nanotechnology 5 protective layers

CE TAP IAT ISO ISO

Better Breathability

Bacterial Filtration Efficiency

Pollen Asbestos Dust Mites Hair Spores

0.07µm

Breathable Fabric

Waterproof Design- Upgrading protection from viruses.

MADE IN TAIWAN

Breathable, secure, hygienic, comfortable, protection, free from stuffiness

190mm

145mm

Schematic diagram of the back of the package (European regulations)

Harmful Particles to human body

									National Certification standard CNS 15719
Pollen 20µm	Dust mite 10µm	Dust 5µm	Charcoal 5µm	PM2.5 2.5µm	Smoke 1µm	Powder 1µm	Incense 0.3µm	Asbestos 0.075µm	

Pass the certification of FFP2 BFE & PFE

Instruction

1. Take out the mask.
2. Place the ear straps on both ears.
3. Press the nose bridge to fit your nose.
4. Adjust the mask in right position.
5. If you feel the mask does not fit well, you can use the slider to correct the two ear straps to make the mask fit the face better.

Caution & Certified

	No smoking		No open flame
	No eating or drinking		No touching the mask
	No touching the mask		CE mark

Notice

- Do not wear masks while exercising.
- Replace the mask once it's torn or twisted.
- Not recommended for those who have asthma, sort of breath, or PMS.

Name: Filtering half mask (Uasmeilve) (TWYC FFP2 MASK)

Applicable object: Adult Model: CMK-16

Specification: 1 Three-dimensional ear hook

2 Mask size 11cm (±0.5cm)*15cm (±0.5cm)

3 Ear strap length 22cm (±0.5cm)

4 Nose piece length 9cm (±0.5cm)

Quantity: 1 piece per package

Executive standard: CNS 1492001+A1:2009 (IEC) 2016435 PPE

Material: PP non-woven fabric, special non-woven fabric support material, polypropylene filter.

Highly protective nanotechnology filter with super soft skin-friendly layer

Lot no: 953000431 Manufacturing date: 2020.06.10

4719871810866

190mm

145mm